



*DURAN® products can be obtained from your specialist laboratory dealer. Order addresses and other details about the items (e. g. technical information) can be found under [www.duran-group.com](http://www.duran-group.com) or in the DURAN® laboratory glassware catalogue.*

*Special safety caps and wide-neck bottles are contained in the brochures "DG Safety Caps" and "DURAN® GLS80 Glass bottles with wide neck".*


## ORDERING MADE EASY:

You would like to equip two bottles with a three port connection system.

You would like to work with two different tube diameters (1.6 mm and 3 mm) and pressure equalization.

Consequently, you need the following individual components:

- Ⓐ Screw cap GL 45 3 port  
(Art. No. 1129751) 1 pack [2 per pack]
- Ⓑ Insert for screw cap GL14  
1.6 mm iØ, (Art. No. 1129815)  
2 packs [1 per pack]
- Ⓒ Screw cap GL 14  
(Art. No. 1129814) 2 packs [2 per pack]
- Ⓓ Insert for screw cap GL14 3.0 mm  
iØ, (Art. No. 1129816) 2 packs [1 per pack]
- Ⓔ Pressure equalization set (complete)  
(Art. No. 1137799) 2 packs [1 per pack]



Order No.	Description	Thread	d mm	h mm	Pack. Qty
11 297 50	Screw cap GL 45, 2-port GL 14	GL 45	54	45	2
11 297 51	Screw cap GL 45, 3-port GL 14	GL 45	54	45	2
11 298 14	Screw cap GL 14 for tube connection	GL 14	20	28	2
11 298 15	Insert for screw cap GL 14, 1.6 mm [1/16 inch] iØ		13	13	1
11 298 16	Insert for screw cap GL 14, 3.0 mm [~1/8 inch] iØ		13	13	1
11 298 17	Insert for screw cap GL 14, 3.2 mm [1/8 inch] iØ		13	13	1
11 298 18	Insert for screw cap GL 14, 6.0 mm [~1/4 inch] iØ		13	13	1
11 562 92	Screw cap GL 14, red	GL 14	20	17	2
11 377 99	Pressure equalization set for 2- and 3-port screw cap (incl. 0.2 µm membrane filter)	GL 14	20	63	1
11 298 19	Spare membrane filter for pressure equalization 0.2 µm				2

Schematic diagram of GL 45 connection system

- ① Screw cap GL 14 (PP)
- ② Silicone sealing ring on insert
- ③ PTFE insert / tubing connector
- ④ Tubing (not supplied)
- ⑤ Port (PP)
- ⑥ O-ring seal

*It's even easier with the DURAN® HPLC laboratory glass bottle, the preassembled complete system*



**DURAN  
WHEATON  
KIMBLE**

Excellence in your hands

DWK Life Sciences GmbH  
Hattenbergstraße 10  
55122 Mainz, Germany  
Tel. +49 (0)6131/1445-4131  
Fax +49 (0)6131/1445-4016  
communications@DWK-LifeSciences.com

[www.DWK-LifeSciences.com](http://www.DWK-LifeSciences.com)

