

MiniMix® 100 P CC

100 mL lab blender

Ref. 011 110 (100-127 V~ / 220-240 V~)

Ref. 011 230 (220-240 V~)



Standards:



Manufactured under:



MiniMix® 100 P CC is designed to blend low setting samples. The adjustable blending power (patented) improves the efficiency of bacterial extraction by adapting to the matrix. MiniMix® 100 P CC is fitted with removable Click & Clean® paddle system (patented) that simplifies cleaning by providing full access to the blending room.

TECHNICAL SPECIFICATIONS

- All stainless steel 304 L
- Multifunction digital screen
- Plain door
- Fully opening door (270°)
- Q-Tight® closing: tight closing for no more leaks from the top of the bag
- Click & Clean® removable paddles (patented)
- Side-by-side paddle stop

- Adjustable blending power (patented)
- Useful volume: 5 - 80 mL
- Variable speed (3 - 6 - 9 - 12 strokes/s)
- Variable blending time (1 s - 1 h or ∞)

- Voltage and frequency
Ref. 011 110 : 100-127 V~ / 220-240 V~ 50/60 Hz
Ref. 011 230 : 220-240 V~ 50/60 Hz
- Fuses
Ref. 011 110 : T3.15A L 250 V / T1.25A L 250 V 5x20 mm
Ref. 011 230 : T1.25A L 250 V 5x20 mm
- Max power: 110 W

- Protection index : 11
- Dimensions (w x d x h): 20 x 37.5 x 20 cm, weight: 11.7 kg
- Box (w x d x h): 50 x 30 x 33 cm, weight: 13.6 kg
- 3-year warranty
- Shock absorbers lifetime warranty
- Availability of spare parts : 10 years

- In compliance with: ISO 7218, ISO 6887-1 and FDA BAM (Bacteriological Analytical Manual)
- Manufactured under CE, RoHS, WEEE
- Designed and made in France

PARIS

Phone: +33 1 34 62 62 61
info@interscience.com

FRANKFURT

Phone: +49 611 7238 7770
sales.germany@interscience.com

BOSTON

Phone: +1 781 937 0007
sales.usa@intersciencelab.com

SHANGHAI

电话: +86 178 2123 6642
sales.china@interscience.cn

SINGAPORE

Phone: +65 6977 7232
sales.asia@interscience.com

TOKYO

Phone: +81 3 6712 9715
sales.japan@interscience.com

interscience

DELIVERED WITH

- 1 power cable
- 1 user's manual
- 1 quick user guide
- 1 free pack of blender bags

OPTIONS / ACCESSORIES

- Ref 016 000 DripTray: Stainless steel security drip tray for MiniMix®

PARIS

Phone: +33 1 34 62 62 61
info@interscience.com

FRANKFURT

Phone: +49 611 7238 7770
sales.germany@interscience.com

BOSTON

Phone: +1 781 937 0007
sales.usa@intersciencelab.com

SHANGHAI

电话: +86 178 2123 6642
sales.china@interscience.cn

SINGAPORE

Phone: +65 6977 7232
sales.asia@interscience.com

TOKYO

Phone: +81 3 6712 9715
sales.japan@interscience.com