



Pipette-/Filter Tips Compatibility Chart

		LLG-proMLP Microliter pipettes							
Art.No.	Volume	6.270 420	6.270 421	6.270 422	6.270 423	6.270 424	6.270 425	6.270 426	6.270 427
		0.2-2 µL	0.5-10 µL	2-20 µL	10-100 µL	20 - 200 µL	100-1000 µL	5-50 mL	1-10 mL

LLG-Pipette Tips Economy 2.0

4.668 775	LLG-Pipette Tips Economy 2.0, 0.1-10 µl, non-sterile, clear, PU=1000	■	■						
4.668 776	LLG-Pipette Tips Economy 2.0, 0,1-10 µl, non-sterile, clear, PU=960	■	■						
4.668 778	LLG-Pipette Tips Economy 2.0, 1-200 µl, non-sterile, clear, PU=960			■	■	■			
4.668 777	LLG-Pipette Tips Economy 2.0, 1-200 µl, non-sterile, yellow, PU=1000			■	■	■			
4.668 779	LLG-Pipette Tips Economy 2.0, 1-200 µl, non-sterile, yellow, PU=960			■	■	■			
4.668 780	LLG-Pipette Tips Economy 2.0, 100-1000 µl, non-sterile, clear, PU=1000						■		
4.668 782	LLG-Pipette Tips Economy 2.0, 100-1000 µl, non-sterile, clear, PU=768						■		
4.679 861	LLG-Pipette Tips Economy 2.0, 100-1000 µl, non-sterile, blue, PU=1000						■		
4.668 783	LLG-Pipette Tips Economy 2.0, 100-1000 µl, non-sterile, blue, PU=768						■		
9.409 097	LLG-Pipette tips Economy 2.0, 100-5000 µl, non-sterile, clear, PU=250							■	
9.409 098	LLG-Pipette tips Economy 2.0, 1000-10.000 µl, non-sterile, clear, PU=100								■

LLG-Filter Tips Typ Low Retention Filter Tips

4.668 784	LLG-Filter Tips 0.1-10 µl Typ Low Retention 2.0, sterile, clear, PU=960	■	■						
4.668 785	LLG-Filter Tips 0.1-10 µl XL Typ Low Retention 2.0, sterile, clear, PU=960	■	■						
4.668 786	LLG-Filter Tips 1-20 µl Typ Low Retention 2.0, sterile, clear, PU=960			■	■	■			
4.668 787	LLG-Filter Tips 1-100 µl Typ Low Retention 2.0, sterile, clear, PU=960				■	■			
4.668 788	LLG-Filter Tips 1-200 µl Typ Low Retention 2.0, sterile, clear, PU=960					■			
4.668 789	LLG-Filter Tips 100-1000 µl Typ Low Retention 2.0, sterile, clear, PU=768						■		

■ compatible

Lab Logistics Group GmbH | Am Hambuch 1 | 53340 Meckenheim | Germany

Fon: +49 (0)2225 - 92 11-0 | Fax: +49 (0)2225 - 92 11-11 | Mail: info@llg-labware.com | Web: www.llg-labware.com